

THE INDUSTRY LEADER IN ALCOHOL BREATH TESTING SINCE 1945

**ALCO-SENSOR IV Total Recall Software**  
**(for the Alco Sensor IV with Memory Instrument)**



**Intoximeters Inc.**

# Alco-Sensor IV Total Recall Software (for the Alco Sensor IV with Memory Instrument)

### GENERAL OVERVIEW

Intoximeters Alco-Sensor IV Total Recall Software combines the precision and accuracy of the Alco-Sensor IV with Memory (ASIV.M) and the processing power of your computer. This software package offers a comprehensive system for managing an alcohol breath test program, storing test data and printing test results.

### ANALYTICAL

Alco-Sensor IV Total Recall Software can interface with supported firmware versions of an ASIV Memory instrument.

### BASIC FUNCTIONALITY

Used with an ASIV.M, the Alco-Sensor IV Total Recall Software collects and stores results from tests run on the ASIV Memory instrument. The software will accept downloads from multiple ASIV.M instruments. The query feature within the database allows the user to search through the stored data, identify records, and generate reports. Reports are generated based on the selection criteria established by the user.

### ALCO-SENSOR IV TOTAL RECALL SOFTWARE FEATURES

#### **Data Storage**

The Alco-Sensor IV Total Recall Software stores tens of thousands of test records, limited only by the disk drive space available on the computer.

#### **Data Collection**

Alco-Sensor IV Total Recall Software requires a serial cable connecting the ASIV.M to the computer for download.

#### **Printout Capability**

Test results and reports can be printed using various printers supported by Windows™ (laser, dot matrix or label).

#### **Communications**

The Alco-Sensor IV Total Recall Software supports a comprehensive on location alcohol breath testing program. Each system can act as either a standalone system supporting a number of instruments or act as part of a larger networked system. Networking multiple Alco-Sensor IV Total Recall Software sites to a centralized database is easily accomplished with IntoxNet, an Intoximeters central database program (purchased separately). Local Alco-Sensor IV Total Recall Software systems can be interfaced with IntoxNet via modem.

#### **Mobility**

Alco-Sensor IV Total Recall Software can be conveniently used wherever your laptop computer goes. If printed results are required, mobile printer options are available.

### MINIMUM CONFIGURATION REQUIRED TO RUN THE ALCO-SENSOR IV TOTAL RECALL SOFTWARE

#### **Central Processing Unit (CPU)**

Pentium class or greater

#### **Random Access Memory (RAM)**

32 mb or greater

#### **Hard Disk Space**

10 mb or greater

#### **Serial Ports (COM Ports)**

1 available port required for cable interface with the Alco-Sensor IV with Memory. A second available port is needed if you choose to use a serial printer.

#### **Parallel Ports (Printer)**

1 required if a parallel printer is used.

#### **Display (Monitor)**

VGA Color or Monochrome Monitor

#### **Input (Keyboard)**

Standard PC keyboard

Supported Magnetic Stripe Readers

#### **Communications (Cable)**

Alco-Sensor IV Total Recall Software comes equipped with an interface cable for the ASIV.M

#### **Communications (Optional)**

9600 or faster baud rate modem.

#### **Operating System**

Microsoft Windows™ 95, 98, ME, NT 4.0 or 2000

## Intoximeters Inc.