

Intoximeters

Experience • Service • Integrity



Alco-Sensor® VxLerator

Downloads breath test results into a Microsoft Excel® compatible format

Intoximeters Alco-Sensor VxLerator is a Windows® application that allows the download of breath test results from the Alco-Sensor VXL into Microsoft Excel® (Intoximeters Alco-Sensor VXL and Microsoft Excel® sold separately).

Alco-Sensor VxLerator will log into your Alco-Sensor VXL and download its stored test results. Click the “View” button, and the data is opened automatically in Microsoft Excel®. You now have the tools of Microsoft Excel® at your disposal to search sort, filter and report on your data.

Alco-Sensor VxLerator benefits (in conjunction with Microsoft Excel®):

- Easy to use
- Allows user to sort multiple breath testing results easily
- Overall organization of your breath testing program
- Makes reporting results easy and seamless
- Provides an organized record-keeping system for past tests

Intoximeters

Data Storage

Maintains a comprehensive record of all testing conducted on your Alco-Sensor VXL. There is no restriction as to how many Alco-Sensor VXL's may be downloaded with a single copy of VXLerator software. Depending on the version of Microsoft Excel, a single spreadsheet can store 1,048,576 (Excel 2007 and beyond) test records.

Ease of Use and Simplicity

One button-click operation turns on the Alco-Sensor VXL and downloads the tests. One button-click opens that data in Excel.

Querying and Statistics

You may use the extensive capabilities native to Microsoft Excel to query and garner statistics from your data.

Minimum Configuration Required to Run the Alco-Sensor Vxlerator System:

Breath Analyzer

Intoximeters Alco-Sensor VXL

PC Operating System

Microsoft Windows 7®, Windows Vista®, or Windows XP® operating system

CPU, Available Hard Disk Space

Requirements do not exceed the minimum system requirements of the underlying Microsoft Windows® operating system.

Microsoft Office Excel Software

Microsoft Office Excel 2010, Excel 2007, Excel 2003, Excel XP, or Excel 2000 (Microsoft Office 9 - 14).

Serial or USB Port

One available USB port or 9-pin serial port is required for connecting your PC to the Alco-Sensor IV. When using a USB port, a supported USB to serial adapter is required (included with package).

Optional Bluetooth Technology

Download available via Bluetooth communication. Bluetooth technology also allows simultaneous multiple instrument download initiation.

Communications Dock and Serial Cable (Included)

Alco-Sensor VXLerator includes a dock and communication cable needed to connect your Alco-Sensor VXL with a standard 9-pin serial port.

Networking

The Alco-Sensor VXLerator Software allows for multi-tiered networking, where data downloaded from multiple workstations is merged at a destination "hub" and into a single spreadsheet, facilitating data consolidation. This networking relies only upon Windows file sharing capability - the availability of a network share (shared directory) common to the workstations and the designated hub.

Microsoft® is a registered trademark of Microsoft® Corp. in the United States and/or other countries.

The names of actual companies and products mentioned herein may be the trademarks of their respective owners.

Alco-Sensor® VXLerator

Downloads breath test results into Microsoft Excel compatible format

